

SCIENTECH®



Setup and Operating Procedures

PN9458P

Thank you for choosing a Scientech Astral calorimeter. Scientech, an ISO 9001 registered company, and our employees are pleased to provide you with a product designed for years of reliable service. Please read this manual completely before using your indicator. This information will enable you to fully utilize the equipment and should be located nearby for reference. The calorimeter is intended to be used only in the manner outlined in this manual. Misuse of the equipment may cause product failure.

The words "calorimeter" and "detector" are synonymous as used in this manual.

Table of Contents:

Calorimeter Operating Parameters:.....	2
Astral™ Calorimeter Specifications:	3
Absorption of HD Absorbing Material:	4
CE Mark Certification:	4
Environmental Requirements:	4
Unpacking and Set Up:	4
Correcting Astral™ Series HD Calorimeters Operating Parameters for Use at Different Wavelengths:	5
Calorimeter Operation Without an Indicator:	5
Cable Requirements:.....	5
Operation of Astral™ Calorimeters with a Digital Volt Meter:.....	6
Operation of Astral™ Calorimeters with an Analog Chart Recorder:	6
Calorimeter Response:	6
Numerical Integration:	7
Initial Voltage Interpolation:.....	7
Peak Voltage Estimate:	7
Calibration of Astral™ Calorimeters Using Electric Substitution Heating:	8
Factory Recalibration:.....	9
Limited Warranty:.....	9
Returned Goods Procedure:	9
Disposal of Electrical and Electronic Equipment:	9
HD Absorption vs. Wavelength:.....	10

CALORIMETER OPERATING PARAMETERS:

Calorimeter 1:

Model No: _____
 Serial No: _____
 Calibration Wavelength _____ nm
 Output Sensitivity (S): _____ V/W
 Time Constant (1/e): _____ sec.
 Calibration Temp: _____ °C
 Sub. Heater Resistance (R_C): _____ ohms
 Sub. Heater Voltage (V_h): _____ volts
 Sub. Heater Wattage (W_h) _____ watts

Calorimeter 2:

Model No: _____
 Serial No: _____
 Calibration Wavelength: _____ nm
 Output Sensitivity (S): _____ V/W
 Time Constant (1/e): _____ sec.
 Calibration Temp: _____ °C
 Sub. Heater Resistance (R_C): _____ ohms
 Sub. Heater Voltage (V_h): _____ volts
 Sub. Heater Wattage (W_h): _____ watts

ASTRAL™ CALORIMETER SPECIFICATIONS:

Model	AC2500	AC25HD	ACX25HD	AC2501	ACX2501	AC25UV	AC2504
Type of Absorber	Surface	Surface	Surface	Volume	Volume	Volume	Volume
Aperture Diameter	25.4 mm	25.4 mm	8 mm	25.4 mm	8 mm	25.4 mm	25.4 mm
Spectral Response	.25-35 μm	.193-12 μm	.4-2 μm	.266-1.2 μm	.4-1.2 μm	.193-.36 μm	.85-4.2 μm
Average Power Maximum	10 W						
Average Power Minimum	1 mW when installed in an Isotherm Enclosure						
Noise Level	10 μW or μJ						
Maximum Power Density	200 W/cm ²	1.5 kW/cm ²	12 kW/cm ²	Note 1	Note 2	Note 3	Note 4
Maximum Peak Power Density	1 MW/cm ²	100 MW/cm ²	800 MW/cm ²	Note 5	8.5 GW/cm ²	Note 6	Note 7
Maximum Single Pulse Energy	10 J						
Maximum Energy Density	Note 8	Note 9	Note 10	Note 11	Note 12	Note 13	Note 14
Precision	< 1 %						
Accuracy	± 3 %						
Response Time	3 sec when connected to a Scientech Indicator in Watts Mode						
Dimensions DxL - inches	3.75 x 2.2	3.75 x 2.2	3.75 x 3.82	3.75 x 2.2	3.75 x 3.82	3.75 x 2.2	3.75 x 2.2
cm	9.53 x 5.6	9.53 x 5.6	9.53 x 9.7	9.53 x 5.6	9.53 x 9.7	9.53 x 5.6	9.53 x 5.6
Weight - pounds/kg	1.5/0.68	1.5/0.68	1.7/0.77	1.5/0.68	1.7/0.77	1.5/0.68	1.5/0.68
Indicator Compatibility	H410, H410D, S310, S310D, D200PC, D200C						

Model	AC5000	AC50HD	ACX50HD	AC5001	ACX5001	AC50UV	AC5004
Type Absorber	Surface	Surface	Surface	Volume	Volume	Volume	Volume
Aperture Diameter	50.8 mm	50.8 mm	16 mm	50.8 mm	16 mm	50.8 mm	50.8 mm
Spectral Response	.25-35 μm	.193-12 μm	.4-2 μm	.266-1.2 μm	.4-1.2 μm	.193-.36 μm	.85-4.2 μm
Average Power Maximum	30 W						
Average Power Minimum	40 mW						
Noise Level	400 μW or μJ						
Maximum Power Density	200 W/cm ²	1.5 kW/cm ²	12 kW/cm ²	Note 1	Note 2	Note 3	Note 4
Maximum Peak Power Density	1 MW/cm ²	100 MW/cm ²	800 MW/cm ²	Note 5	8.5 GW/cm ²	Note 6	Note 7
Maximum Single Pulse Energy	30 J						
Maximum Energy Density	Note 8	Note 9	Note 10	Note 11	Note 12	Note 13	Note 14
Precision	< 1 %						
Accuracy	± 3 %						
Response Time	3 sec when connected to a Scientech Indicator in Watts Mode						
Dimensions DxL - inches	4.75 x 2.3	4.75 x 2.3	4.75 x 3.92	4.75 x 2.3	4.75 x 3.92	4.75 x 2.3	4.75 x 2.3
cm	12.07 x 5.8	12.07 x 5.8	12.07 x 9.96	12.07 x 5.8	12.07 x 9.96	12.07 x 5.8	12.07 x 5.8
Weight pounds/kgs	2.9/1.3	2.9/1.3	3.1/1.4	2.9/1.3	3.1/1.4	2.9/1.3	2.9/1.3
Indicator Compatibility	H410, H410D, S310, S310D, D200PC, D200C						

Note 1: AC2501, AC5001 30W/cm² @ 1064nm, 23W/cm² @ 532nm, 8.5W/cm² @ 355nm, 175mW/cm² @ 266nm

Note 2: ACX2501, ACX5001 Note 1 specs x 8 for 400nm to 1.2μm

Note 3: AC25UV, AC50UV 50W/cm² @ 355nm

Note 4: AC2504, AC5004 35W/cm² @ 1064nm

Note 5: AC2501, AC5001 100GW/cm² @ 1064nm, 78GW/cm² @ 532nm, 29GW/cm² @ 355nm, 580MW/cm² @ 266nm

Note 6: AC25UV, AC50UV Repetitive pulses: 101MW/cm² @ 355nm

Single pulses: 3.5GW/cm² @ 355nm

Note 7: AC2504, AC5004 125GW/cm² @ 1064nm

Note 8: AC2500, AC5000 Max J/cm² = 1,000 x (pulse width)^{1/2} to a maximum of 200J/cm².

Note 9: AC25HD, AC50HD Max J/cm² = 4,500 x (pulse width)^{1/2} to a maximum of 14J/cm².

Note 10: ACX25HD, ACX50HD Max J/cm² = 36,000 x (pulse width)^{1/2} to a maximum of 42.5J/cm².

Note 11: AC2501, AC5001 Repetitive pulses: 4.1J/cm²@1064nm, 3.2J/cm²@532nm, 1.2J/cm²@355nm, 24mJ/cm²@266nm

Single pulses: 8J/cm²@1064nm, 6.2J/cm²@532nm, 2.3J/cm²@355nm, 46mJ/cm²@266nm

Note 12: ACX2501, ACX5001 Note 11 specs x 8 for 400nm to 1.2μm

Note 13: AC25UV, AC50UV Repetitive pulses: 1.1J/cm² @ 355nm

Single pulses: 40J/cm² @ 355nm

Note 14: AC2504, AC5004 Repetitive pulses: 4.8J/cm² @ 1064nm

Single pulses: 10J/cm² @ 1064nm

ABSORPTION OF HD ABSORBING MATERIAL:

Warning: You must exercise caution when using HD detectors. They exhibit spectral reflection of between 7% and 18%, of the input power, back out of the aperture. Please refer to Figure 1 to determine the reflectance for the wavelength you are measuring. These detectors should be treated as a partial mirror or any other type of reflective optic and the appropriate caution level observed, especially at the CO₂ wavelength.




Figure 1

Note: HD calorimeters might show “beam” marks on the absorbing surface. These marks are characteristic of the material and do not affect the performance of the calorimeter. Do not rub off or remove the marks. Polishing or cleaning the absorbing surface might change the performance of the calorimeter.

Please see the table at the end of this manual for specific absorption vs. wavelength information.

CE MARK CERTIFICATION:

All of the calorimeters listed in this manual have been certified for the European CE mark.

ENVIRONMENTAL REQUIREMENTS:

This product is intended for indoor use at altitudes up to 2000 meters, Pollution Degree 2 in accordance with IEC 664 and transient overvoltages according to Installation Categories (Overvoltage Categories) II. Note that each of the above detectors will not pass the IEC 801 Publication, Part 3, Radiated Electromagnetic Field Requirements. The system, meter and detector, is designed to measure radiation within the test's radiation band. The detectors were held outside the radiated electromagnetic field during this test. It is up to the user to be aware of RF fields present during measurements and their effects if any on those measurements.

UNPACKING AND SET UP:

The calorimeter and accessories are shipped in custom packing materials. All packing materials should be saved for future damage free shipments.

A ½” diameter mounting post is included. Screw the post into one of the mounting holes in the body of the calorimeter. Mount the post to your optical bench or working surface. An optional mounting base, Scientech Model 301-019, is available for this purpose.

A 3 meter mini DIN interconnect cable is also included. If you are using a Scientech indicator, connect the calorimeter to the indicator with the interconnect cable. Follow the detailed set up instructions that are in the indicator's instruction manual. If you are not using a Scientech indicator operating requirements are contained in this manual.

Note: Astral calorimeters are sensitive to all types of thermal input. Due to the handling of the calorimeter during setup and possible environmental temperature differences, thermal gradients may exist in the calorimeter. Allow the calorimeter to sit undisturbed for several minutes to reach thermal equilibrium, before using.

Note: When using a 25mm Astral calorimeter for measuring average power levels below 30mW and single pulse energy levels below 30mJ, a Scientech Model 36-0203A, Isoperibol Enclosure, is highly recommended. The isoperibol enclosure should not be used at average power levels above 30mW, and single pulse energy levels above 100mJ because heat build up will occur.

CORRECTING ASTRAL™ SERIES HD CALORIMETERS OPERATING PARAMETERS FOR USE AT DIFFERENT WAVELENGTHS:

Scientech Astral calorimeters in general have a flat response to all wavelengths within their specified spectral response. HD and HDX calorimeters are an exception to that rule and are calibrated at a specific wavelength by adjusting the calorimeter's gain circuitry for that wavelength. The calibration wavelength is recorded in the Operating Parameters section at the front of the manual and on the detector's serial tag. When a HD or HDX calorimeter is used at a wavelength other than the calibration wavelength, the calorimeter's output sensitivity can be adjusted to compensate for the absorption rate at the new wavelength as follows:

1. Find the absorption rate from the chart at the end of this manual for the calibration wavelength of your calorimeter.
2. Find the absorption rate for the wavelength where you will be working.
3. Determine the new output sensitivity using the following formula:

$$\frac{\text{absorption rate of calibration wavelength}}{\text{absorption rate of new wavelength}} \times \text{output sensitivity from serial tag} = \text{output sensitivity for new wavelength}$$

NOTE: Due to variability in the manufacturing process the absorption characteristics of the HD model calorimeters can vary in the UV region (190 to 400nm). Scientech recommends optical calibration at 266nm if the detector is to be used in the UV region instead of relying on an absorption rate in the chart at the end of this manual.

CALORIMETER OPERATION WITHOUT AN INDICATOR:

Cable Requirements:

Astral calorimeters are powered up by the indicators. To use an Astral calorimeter without a Scientech indicator, but with a volt meter or chart recorder, you must apply +/-8VDC to the mini DIN connector as shown in Figure 2. The voltage output of the calorimeter, from pin 8, should be connected to the positive side of the DVM or chart recorder. All 3 of the grounds should be tied together at the negative side. Pins 2 and 3 are not used.




Figure 2

Operation of Astral™ Calorimeters with a Digital Volt Meter:

The calorimeters may be used with any digital volt meter (DVM) capable of reading 5 volts full scale.

Connect the output of the calorimeter to the DVM.

Select the DC volts mode.

Direct the laser beam on to the absorbing surface of the calorimeter.

When the display of the DVM has stabilized (about 2 minutes), calculate the laser power using the formula:

$$W = V/S$$

where:

W = Laser power in watts

V = Voltage reading of the DVM in volts

S = Sensitivity of the calorimeter from page 2.

Operation of Astral™ Calorimeters with an Analog Chart Recorder:

Calorimeter Response:

The response of a calorimeter to a single pulse input as displayed by a chart recorder appears below.




Figure 3

The output voltage from a chart recorder can be converted to wattage at any time by:

$$W = V/S, W_i = V_i/S$$

V = Chart recorder voltage level in mV

S = Calorimeter sensitivity in mV/W

The total energy (E) in the pulse can be found by integrating the instantaneous wattage over time:

$$E = \int_0^{\infty} W(t) dt$$

The following methods may be used to compute the total integrated energy:

Numerical Integration:

Finding the area under the curve in figure 7 is the equivalent procedure for determining pulse energy. Choose an appropriate time interval, dt , and perform the summation:

$$E = \sum_{i=1}^N W_i \times dt = (dt/S) \sum_{i=1}^N V_i$$

The error caused by this procedure is:

$$dE = (dt/S) \sum_{i=1}^N dV_i$$

The error, in theory, is only dependent upon the value of $\sum dV_i$, that is the cumulative random error of V_i . This number should approach zero if data is carefully taken. The accuracy is also increased if the time interval, dt , is minimized. Numerical integration can yield accurate results, but is a tedious task.

Initial Voltage Interpolation:

A method used to eliminate the tedious numerical integration task is to project the thermal decay envelope on to the voltage axis, determine the $1/e$ decay time constant T , and estimate the total energy value (E):

$$E = (V_o/S) \times T$$

The change from thermal absorption to thermal transport phenomena near the peak causes difficulty in accurately projecting the envelope on to the voltage axis introducing an error, dV_o . Further, the determination of the time constant T , introduces another error, dT . The total error is the sum of the two errors.

$$dE = (V_o/S)dT + (T/S)dV_o$$

The difficulty in eliminating the potential error makes this method typically less accurate than numerical integration, but much faster in application.

Peak Voltage Estimate:

The peak voltage method requires using an independent determination of total energy and referencing it back to the peak voltage value, V_p .

For a given pulse, use the numerical integration method to obtain E . Note the peak voltage, V_p . Compute the value, F

$$F = E/V_p$$

For the next pulse compute the total energy: $E = F \times V_p$

The error in using this method yields: $dE = FdV_p + V_p dF$

The accuracy of this measurement depends upon the error in the original calibration, dF , and the error in the peak voltage dV_p . A careful numerical integration yields a value for dF near zero. The value of dV_p can be minimized by maintaining the geometry of the system (i.e. beam intensity, beam profile, wavelength and environment) during operation to be the same as during calibration. Under controlled circumstances, the peak method accuracy usually falls between the numerical integration and initial voltage interpolation methods.

CALIBRATION OF ASTRAL™ CALORIMETERS USING ELECTRIC SUBSTITUTION HEATING:

For Astral calorimeters the electric substitution heating option must be ordered and installed at the factory when the calorimeter is purchased. It can not be retrofitted to a calorimeter at a later time. To calibrate using electric substitution heating proceed as follows.

- A. Remove the screws holding the calorimeter's ID tag and remove the plate to expose the circuit board.

Calorimeter Circuit Board




Figure 4

- B. Connect a DVM to the test points labeled SUB and HTR on the calorimeter circuit board.
- C. Measure the resistance of the substitution heater making sure to subtract the resistance of the patch cables from the total resistance measurement. Compare this resistance to R_c in the calibration data in the front of the manual. The two should agree within 2%. If not contact Scientech.
- D. Remove the DVM. Connect a power supply to the SUB and HTR test points and connect the DVM to monitor the power supply.
- E. Connect a second DVM to the output and ground of the Astral calorimeter.
- F. Apply V_h volts, stated in the calibration data you received with the calorimeter, to the substitution heater.
- G. Verify that the output voltage is one of the following values depending on the type of calorimeter being tested:
- 0.750 volts for 25 mm calorimeters
 - 0.333 volts for 50 mm calorimeters
- H. If needed, adjust the calibration trim pot, R4 on the calorimeter circuit board, until the voltage measured in step G is from the calibration data, is displayed by the indicator.

FACTORY RECALIBRATION:

Scientech recommends that a complete calibration be performed annually to verify system accuracy. Please contact our Product Service Department at (800)525-0522 or (303)444-1361 or Fax (303)444-9229 or email inst@scientech-inc.com to arrange for a NIST traceable, factory calibration.

LIMITED WARRANTY:

All Scientech Laser Power and Energy Measurement Systems are warranted against defects in materials and workmanship for two (2) years from date of delivery. During the warranty period, Scientech will repair, or at its option replace at no charge, components that prove to be defective. The equipment must be returned, shipping prepaid, to Scientech's product service facility. This limited warranty does not apply if the equipment is damaged by accident or misuse or as a result of service or modification by other than a Scientech service facility. The foregoing warranty is in lieu of all other warranties expressed or implied including but not limited to any implied warranty of merchantability, fitness, or adequacy for any special incidental or consequential damages whether in contract, tort, or otherwise.

RETURNED GOODS PROCEDURE:

Should it become necessary to return any item to Scientech for any reason, please contact our Product Service Department at (800)525-0522 or (303)444-1361 or Fax (303)444-9229 or email inst@scientech-inc.com. When you call, please be ready to provide model number, serial number, and a description of the problem. Frequently we can provide self-help information which will eliminate the need for returning the unit(s).


If equipment return is required, please pack the items in the original box and packing material. As an alternate, place the equipment in a snug-fitting box, and then pack that box in a larger box with at least four inches of packing material. Scientech does not assume responsibility for under packed items. Please include the name and phone number of the person we should contact regarding repair questions. Normally, products are repaired and shipped within 5 working days after their arrival at the product service facility. This is an average time and could vary depending on the workload.

Shipping address: Scientech, Inc.
 Product Service Department
 5649 Arapahoe Ave.
 Boulder, Colorado 80303 U.S.A.

DISPOSAL OF ELECTRICAL AND ELECTRONIC EQUIPMENT:

Scientech, Inc. recommends the following for disposal of electrical and electronic equipment:

1. The best option is to reuse the equipment in its entirety.
2. Where the equipment can not be reused in its entirety, priority should be given to reuse of its subassemblies and components.
3. Where reuse is not appropriate, electrical and electronic equipment, including batteries, should be recycled according to local ordinances.
4. Waste electrical and electronic equipment should never be mixed with municipal waste.



HD ABSORPTION VS. WAVELENGTH:

Wavelength μm	Absorption %
0.200	93.66
0.210	93.78
0.220	92.18
0.230	90.44
0.240	88.79
0.250	90.07
0.255	90.42
0.260	90.41
0.265	90.56
0.270	90.71
0.275	91.13
0.280	91.63
0.285	91.85
0.290	92.01
0.295	92.73
0.300	93.02
0.305	93.04
0.310	93.15
0.315	93.53
0.320	93.45
0.325	93.82
0.330	93.83
0.335	93.96
0.340	94.07
0.345	94.14
0.350	94.15
0.355	94.19
0.360	94.31
0.365	94.30
0.370	94.44
0.375	94.48
0.380	94.46
0.385	94.53
0.390	94.53
0.395	94.57
0.400	94.62
0.405	94.63
0.410	94.65
0.415	94.72
0.420	94.70
0.425	94.78
0.430	94.74
0.435	94.80
0.440	94.76
0.445	94.80
0.450	94.79
0.455	94.79
0.460	94.79
0.465	94.76
0.470	94.78
0.475	94.75

Wavelength μm	Absorption %
0.480	94.71
0.485	94.75
0.490	94.80
0.495	94.76
0.500	94.73
0.505	94.70
0.510	94.67
0.515	94.68
0.520	94.70
0.525	94.66
0.530	94.66
0.535	94.70
0.540	94.75
0.545	94.73
0.550	94.66
0.555	94.63
0.560	94.64
0.565	94.64
0.570	94.62
0.575	94.56
0.580	94.59
0.585	94.56
0.590	94.56
0.595	94.61
0.600	94.55
0.605	94.51
0.610	94.50
0.615	94.45
0.620	94.38
0.625	94.44
0.630	94.42
0.635	94.44
0.640	94.48
0.645	94.50
0.650	94.46
0.655	94.49
0.660	94.39
0.665	94.43
0.670	94.41
0.675	94.44
0.680	94.36
0.685	94.37
0.690	94.36
0.695	94.19
0.700	94.37
0.705	94.25
0.710	94.20
0.715	94.16
0.720	94.08
0.725	94.11
0.730	94.21

Wavelength μm	Absorption %
0.735	93.99
0.740	94.21
0.745	93.98
0.750	93.98
0.755	94.08
0.760	94.04
0.765	93.89
0.770	93.94
0.775	94.01
0.780	93.91
0.785	94.03
0.790	93.96
0.795	93.89
0.800	93.64
0.805	93.65
0.810	93.67
0.815	93.53
0.820	94.21
0.825	93.61
0.830	93.81
0.835	93.83
0.840	93.83
0.845	93.93
0.850	93.84
0.855	93.31
0.860	93.65
0.865	93.90
0.870	93.84
0.875	93.92
0.880	93.81
0.885	93.85
0.890	93.54
0.895	93.88
0.900	93.60
0.905	93.70
0.910	93.62
0.915	93.68
0.920	93.66
0.925	93.64
0.930	93.60
0.935	93.67
0.940	93.59
0.945	93.60
0.950	93.50
0.955	93.49
0.960	93.51
0.965	93.45
0.970	93.42
0.975	93.40
0.980	93.35
0.985	93.37

Wavelength μm	Absorption %
0.990	93.34
0.995	93.35
1.000	93.34
1.005	93.33
1.010	93.33
1.015	93.33
1.020	93.29
1.025	93.28
1.030	93.28
1.035	93.25
1.040	93.23
1.045	93.17
1.050	93.14
1.055	93.12
1.060	93.12
1.065	93.14
1.070	93.17
1.075	93.16
1.080	93.13
1.085	93.12
1.090	93.09
1.095	93.06
1.100	93.03
1.105	93.01
1.110	93.00
1.115	92.83
1.120	92.76
1.125	92.67
1.130	92.60
1.135	92.57
1.140	92.56
1.145	92.51
1.150	92.47
1.155	92.43
1.160	92.42
1.165	92.37
1.170	92.32
1.175	92.30
1.180	92.24
1.185	92.23
1.190	92.21
1.195	92.16
1.200	92.08
1.205	92.02
1.210	91.94
1.215	91.91
1.220	91.92
1.225	91.88
1.230	91.87
1.235	91.82
1.240	91.81

Wavelength μm	Absorption %
1.245	91.78
1.250	91.79
1.255	91.76
1.260	91.75
1.265	91.76
1.270	91.79
1.275	91.75
1.280	91.74
1.285	91.74
1.290	91.74
1.295	91.71
1.300	91.72
1.305	91.70
1.310	91.70
1.315	91.73
1.320	91.78
1.325	91.78
1.330	91.79
1.335	91.76
1.340	91.77
1.345	91.82
1.350	91.94
1.355	92.17
1.360	92.17
1.365	92.17
1.370	92.15
1.375	92.12
1.380	92.12
1.385	92.12
1.390	92.10
1.395	92.04
1.400	91.99
1.405	92.00
1.410	92.02
1.415	91.98
1.420	91.95
1.425	91.95
1.430	91.96
1.435	91.94
1.440	91.91
1.445	91.86
1.450	91.77
1.455	91.77
1.460	91.80
1.465	91.79
1.470	91.79
1.475	91.77
1.480	91.75
1.485	91.73
1.490	91.69
1.495	91.69
1.500	91.68
1.505	91.71

Wavelength μm	Absorption %
1.510	91.70
1.515	91.66
1.520	91.58
1.525	91.59
1.530	91.59
1.535	91.58
1.540	91.57
1.545	91.53
1.550	91.52
1.555	91.50
1.560	91.49
1.565	91.51
1.570	91.50
1.575	91.49
1.580	91.49
1.585	91.50
1.590	91.47
1.595	91.45
1.600	91.44
1.605	91.44
1.610	91.43
1.615	91.36
1.620	91.32
1.625	91.33
1.630	91.34
1.635	91.34
1.640	91.31
1.645	91.29
1.650	91.24
1.655	91.26
1.660	91.30
1.665	91.24
1.670	91.24
1.675	91.20
1.680	91.15
1.685	91.16
1.690	91.18
1.695	91.13
1.700	91.09
1.705	91.04
1.710	91.03
1.715	90.97
1.720	90.99
1.725	90.96
1.730	90.92
1.735	90.89
1.740	90.92
1.745	90.92
1.750	90.89
1.755	90.89
1.760	90.89
1.765	90.86
1.770	90.87

Wavelength μm	Absorption %
1.775	90.84
1.780	90.86
1.785	90.83
1.790	90.78
1.795	90.77
1.800	90.80
1.805	90.79
1.810	90.77
1.815	90.76
1.820	90.75
1.825	90.76
1.830	90.81
1.835	90.71
1.840	90.74
1.845	90.71
1.850	90.77
1.855	90.71
1.860	90.74
1.865	90.70
1.870	90.73
1.875	90.67
1.880	90.67
1.885	90.64
1.890	90.64
1.895	90.61
1.900	90.59
1.905	90.53
1.910	90.48
1.915	90.63
1.920	90.45
1.925	90.48
1.930	90.48
1.935	90.48
1.940	90.52
1.945	90.45
1.950	90.41
1.955	90.47
1.960	90.50
1.965	90.56
1.970	90.47
1.975	90.49
1.980	90.46
1.985	90.36
1.990	90.39
1.995	90.36
2.000	90.39
2.001	90.68
2.004	90.79
2.005	90.49
2.007	90.53
2.010	90.56
2.013	90.54
2.015	90.52

Wavelength μm	Absorption %
2.016	90.87
2.019	91.03
2.020	90.69
2.022	90.68
2.025	90.67
2.026	90.99
2.029	90.92
2.030	90.60
2.032	90.61
2.035	90.63
2.038	90.66
2.040	90.69
2.041	90.70
2.045	90.70
2.048	90.69
2.050	90.68
2.051	91.03
2.054	91.07
2.055	90.70
2.058	90.67
2.060	90.63
2.061	91.02
2.064	91.03
2.065	90.64
2.068	90.64
2.070	90.63
2.071	91.01
2.074	90.96
2.075	90.55
2.077	90.55
2.080	90.54
2.081	90.91
2.084	90.76
2.085	90.43
2.088	90.47
2.090	90.51
2.091	90.87
2.094	90.90
2.095	90.53
2.098	90.57
2.100	90.61
2.101	90.99
2.104	90.90
2.105	90.50
2.108	90.45
2.110	90.40
2.111	90.41
2.115	90.43
2.118	90.80
2.120	89.36
2.125	90.44
2.122	90.44
2.125	90.60

Wavelength μm	Absorption %
2.129	90.76
2.130	90.41
2.132	90.05
2.135	90.41
2.136	90.43
2.139	90.84
2.140	90.48
2.143	90.46
2.145	91.65
2.146	90.04
2.150	90.45
2.153	90.44
2.155	90.42
2.157	90.39
2.160	90.36
2.161	90.69
2.164	90.76
2.165	90.43
2.168	90.44
2.170	90.44
2.171	90.41
2.175	90.38
2.179	90.38
2.180	90.37
2.182	90.40
2.185	90.44
2.186	90.42
2.190	90.41
2.193	90.39
2.195	90.37
2.197	90.38
2.200	90.39
2.201	90.34
2.205	90.29
2.208	90.33
2.210	90.36
2.212	90.33
2.215	90.29
2.216	90.27
2.220	90.26
2.224	90.21
2.225	90.16
2.227	90.20
2.230	90.24
2.231	90.24
2.235	90.25
2.239	90.22
2.240	90.19
2.243	90.14
2.245	90.08
2.247	90.17
2.250	90.25
2.251	90.21

Wavelength μm	Absorption %
2.255	90.17
2.258	90.14
2.260	90.11
2.262	90.12
2.265	90.13
2.270	90.24
2.274	90.20
2.275	90.17
2.278	90.17
2.280	90.17
2.282	90.14
2.285	90.11
2.286	90.15
2.290	90.19
2.294	90.11
2.295	90.03
2.298	90.10
2.300	90.18
2.303	90.14
2.305	90.10
2.307	90.08
2.310	90.07
2.311	90.09
2.315	90.10
2.319	90.07
2.320	90.04
2.323	90.03
2.325	90.02
2.327	90.00
2.330	89.98
2.332	89.97
2.335	89.95
2.336	89.94
2.340	89.93
2.344	90.08
2.345	90.23
2.348	90.08
2.350	89.93
2.353	89.96
2.355	89.99
2.357	89.98
2.360	89.98
2.361	89.88
2.365	89.79
2.366	89.82
2.370	89.85
2.374	89.98
2.375	90.10
2.379	89.97
2.380	89.84
2.383	89.82
2.385	89.79
2.387	89.86

Wavelength μm	Absorption %
2.390	89.92
2.392	89.92
2.395	89.91
2.396	89.72
2.400	89.53
2.401	89.57
2.405	89.61
2.410	89.84
2.414	89.79
2.415	89.73
2.419	89.85
2.420	89.96
2.423	89.88
2.425	89.79
2.428	89.91
2.430	90.03
2.432	89.89
2.435	89.75
2.437	89.68
2.440	89.61
2.441	89.78
2.445	89.95
2.446	89.90
2.450	89.86
2.451	89.73
2.455	89.59
2.460	89.75
2.465	89.83
2.469	89.81
2.470	89.80
2.474	89.66
2.475	89.51
2.479	89.72
2.480	89.92
2.483	89.79
2.485	89.66
2.488	89.82
2.490	89.98
2.493	89.79
2.495	89.59
2.498	89.49
2.500	89.38
2.503	89.38
2.507	89.38
2.512	89.35
2.517	89.37
2.522	89.39
2.527	89.41
2.532	89.43
2.537	89.36
2.542	89.34
2.547	89.27
2.552	89.24

Wavelength μm	Absorption %
2.557	89.30
2.562	89.25
2.567	89.23
2.572	89.29
2.577	89.25
2.582	89.28
2.588	89.26
2.593	89.20
2.598	89.18
2.603	89.18
2.608	89.13
2.614	89.11
2.619	89.12
2.624	89.16
2.630	89.10
2.635	89.10
2.640	89.13
2.646	89.20
2.651	89.21
2.656	89.11
2.662	89.06
2.667	89.03
2.673	88.98
2.678	89.03
2.684	89.06
2.690	89.10
2.695	89.05
2.701	89.09
2.706	89.10
2.712	89.05
2.718	89.06
2.723	89.05
2.729	89.09
2.735	89.09
2.741	89.13
2.747	89.19
2.752	89.19
2.758	89.16
2.764	89.20
2.770	89.23
2.776	89.23
2.782	89.24
2.788	89.30
2.794	89.32
2.800	89.27
2.806	89.27
2.812	89.26
2.818	89.25
2.824	89.27
2.830	89.26
2.837	89.24
2.843	89.29
2.849	89.31

Wavelength μm	Absorption %
2.855	89.26
2.862	89.24
2.868	89.25
2.874	89.23
2.881	89.24
2.887	89.26
2.894	89.24
2.900	89.19
2.907	89.22
2.913	89.20
2.920	89.22
2.926	89.24
2.933	89.20
2.940	89.21
2.946	89.17
2.953	89.17
2.960	89.17
2.966	89.17
2.973	89.15
2.980	89.12
2.987	89.14
2.994	89.17
3.001	89.14
3.008	89.11
3.015	89.09
3.022	89.08
3.029	89.08
3.036	89.06
3.043	89.07
3.050	89.06
3.057	89.04
3.065	89.03
3.072	89.02
3.079	89.01
3.087	89.00
3.094	89.01
3.101	89.00
3.109	88.99
3.116	88.99
3.124	88.95
3.131	88.93
3.139	88.97
3.146	88.96
3.154	88.92
3.162	88.89
3.170	88.88
3.177	88.86
3.185	88.85
3.193	88.86
3.201	88.85
3.209	88.83
3.217	88.82
3.225	88.80

Wavelength μm	Absorption %
3.233	88.79
3.241	88.75
3.249	88.75
3.257	88.76
3.265	88.75
3.274	88.74
3.282	88.73
3.290	88.72
3.299	88.71
3.307	88.69
3.315	88.68
3.324	88.67
3.333	88.65
3.341	88.65
3.350	88.63
3.358	88.62
3.367	88.61
3.376	88.59
3.385	88.57
3.394	88.55
3.402	88.51
3.411	88.35
3.420	88.05
3.429	87.89
3.439	88.03
3.448	88.18
3.457	88.23
3.466	88.26
3.475	88.29
3.485	88.30
3.494	88.28
3.504	88.11
3.513	88.01
3.523	88.11
3.532	88.19
3.542	88.23
3.552	88.26
3.561	88.24
3.571	88.21
3.581	88.22
3.591	88.22
3.601	88.22
3.611	88.21
3.621	88.21
3.631	88.20
3.641	88.20
3.652	88.19
3.662	88.18
3.672	88.17
3.683	88.15
3.693	88.14
3.704	88.13
3.714	88.11

Wavelength μm	Absorption %
3.725	88.11
3.736	88.10
3.747	88.09
3.758	88.07
3.768	88.07
3.779	88.07
3.790	88.06
3.802	88.05
3.813	88.04
3.824	88.02
3.835	88.03
3.847	88.02
3.858	88.00
3.870	87.99
3.881	87.99
3.893	87.98
3.905	87.97
3.916	87.96
3.928	87.96
3.940	87.96
3.952	87.95
3.964	87.94
3.977	87.93
3.989	87.90
4.001	87.90
4.013	87.90
4.026	87.88
4.038	87.89
4.051	87.89
4.064	87.90
4.077	87.88
4.089	87.85
4.102	87.85
4.115	87.85
4.129	87.83
4.142	87.82
4.155	87.82
4.168	87.82
4.182	87.80
4.195	87.77
4.209	87.70
4.223	87.60
4.236	87.56
4.250	87.60
4.264	87.60
4.278	87.57
4.293	87.58
4.307	87.61
4.321	87.64
4.336	87.65
4.350	87.66
4.365	87.67
4.380	87.67

Wavelength μm	Absorption %
4.394	87.66
4.409	87.64
4.424	87.63
4.440	87.62
4.455	87.61
4.470	87.61
4.486	87.59
4.501	87.58
4.517	87.56
4.533	87.55
4.549	87.55
4.565	87.54
4.581	87.51
4.597	87.49
4.613	87.50
4.630	87.49
4.646	87.48
4.663	87.46
4.680	87.45
4.697	87.44
4.714	87.43
4.731	87.41
4.749	87.40
4.766	87.39
4.784	87.38
4.801	87.37
4.819	87.36
4.837	87.35
4.855	87.34
4.873	87.32
4.892	87.31
4.910	87.30
4.929	87.29
4.948	87.27
4.967	87.25
4.986	87.24
5.005	87.22
5.025	87.20
5.044	87.20
5.064	87.18
5.084	87.16
5.104	87.14
5.124	87.13
5.144	87.10
5.165	87.10
5.185	87.09
5.206	87.08
5.227	87.07
5.248	87.05
5.270	87.03
5.291	87.01
5.313	87.00
5.335	86.98

Wavelength μm	Absorption %	Wavelength μm	Absorption %	Wavelength μm	Absorption %	Wavelength μm	Absorption %
5.357	86.95	6.682	86.05	8.879	84.32	13.09	80.67
5.379	86.95	6.717	86.06	8.940	84.27	13.23	80.54
5.401	86.94	6.752	86.07	9.002	84.22	13.36	80.36
5.424	86.90	6.787	86.00	9.065	84.17	13.50	80.24
5.447	86.89	6.823	85.91	9.129	84.11	13.65	80.08
5.470	86.86	6.859	85.86	9.194	84.06	13.79	79.83
5.493	86.86	6.895	85.91	9.260	84.00	13.94	79.67
5.516	86.86	6.932	85.93	9.326	83.95	14.09	79.57
5.540	86.83	6.970	85.92	9.394	83.90	14.25	79.42
5.564	86.80	7.007	85.91	9.462	83.83	14.40	79.29
5.588	86.78	7.045	85.87	9.532	83.78	14.57	79.14
5.612	86.79	7.084	85.86	9.603	83.72	14.73	78.98
5.636	86.77	7.123	85.84	9.674	83.64	14.90	78.81
5.661	86.76	7.162	85.79	9.747	83.55	15.07	78.66
5.686	86.76	7.202	85.77	9.821	83.46	15.25	78.54
5.711	86.74	7.242	85.74	9.896	83.38	15.43	78.39
5.736	86.68	7.283	85.70	9.972	83.32	15.62	78.23
5.762	86.59	7.324	85.66	10.05	83.26	15.81	78.05
5.787	86.62	7.366	85.63	10.13	83.19	16.00	77.89
5.813	86.66	7.408	85.61	10.21	83.13	16.20	77.73
5.839	86.66	7.450	85.56	10.29	83.06	16.41	77.54
5.866	86.64	7.493	85.52	10.37	82.99	16.62	77.36
5.893	86.61	7.537	85.50	10.45	82.90	16.84	77.17
5.919	86.63	7.581	85.48	10.54	82.82	17.06	76.95
5.947	86.63	7.626	85.46	10.60	82.79	17.29	76.77
5.974	86.64	7.671	85.42	10.63	82.76	17.52	76.57
6.002	86.67	7.716	85.36	10.71	82.69	17.76	76.38
6.030	86.66	7.763	85.32	10.80	82.61	18.01	76.15
6.058	86.60	7.809	85.28	10.89	82.51	18.26	75.87
6.086	86.62	7.857	85.24	10.99	82.42	18.52	75.62
6.115	86.64	7.905	85.18	11.08	82.35	18.79	75.36
6.144	86.63	7.953	85.13	11.18	82.29	19.06	75.07
6.173	86.62	8.002	85.09	11.27	82.21	19.35	74.81
6.203	86.62	8.052	85.02	11.37	82.13	19.64	74.56
6.232	86.61	8.102	84.97	11.47	82.06	19.94	74.25
6.263	86.58	8.153	84.94	11.58	81.94	20.26	73.99
6.293	86.54	8.205	84.91	11.68	81.86	20.58	73.73
6.324	86.47	8.257	84.86	11.79	81.79	20.91	73.41
6.355	86.39	8.310	84.82	11.89	81.69	21.25	73.14
6.386	86.34	8.364	84.77	12.00	81.60	21.61	72.90
6.418	86.26	8.418	84.73	12.12	81.50	21.97	72.65
6.450	86.23	8.473	84.69	12.23	81.40	22.35	72.37
6.482	86.15	8.529	84.60	12.35	81.29	22.74	71.98
6.514	86.12	8.585	84.50	12.47	81.21	23.15	71.60
6.547	86.10	8.642	84.47	12.59	81.10	23.57	71.31
6.580	86.05	8.700	84.45	12.71	80.98	24.01	70.85
6.614	86.03	8.759	84.41	12.84	80.85	24.46	70.50
6.648	86.03	8.819	84.37	12.96	80.74		

Note: Due to variability in the manufacturing process the absorption in the UV region varies.
Calibration at 266nm is recommended if the detector is to be used in the UV region.